


NEPHROSTOMY KIT



Set components:

- Pigtail type catheter
- guidewire "J"
- two-part needle or three-part needle
- dilators
- dilator with sleeve to peel away
- fixing collar
- syringe 10 ml LL
- scalpel
- band clip

The set is intended for urine diversion through nephrostomy opening in case of blockade of natural urine outflow from the kidney.

Reference:	Catheter size [F]	Catheter length [cm]	Dilator with sleeve to peel away size [F]	Needle type	Needle size/ length	Guidewire diameter [inch]
ZNEF5F	5	45	5, 6	Three-part needle	20 Gx19 cm	0.035"
ZNEF6F	6		7, 8			
ZNEF8F	8		7, 10			
ZNEF9F	9		7, 10	Two-part needle	18 Gx20 cm	0.038"
ZNEF10F	10		7, 10, 13			
ZNEF11F	11		7, 10, 13			
ZNEF12F	12		7, 10, 13			
ZNEF14F	14		10, 13, 15			

When ordering use Ref. Number, e.g.: ZNEF14F.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

