




# Urology

Marflow AG  
Soodstrasse 57  
8134 Adliswil/Zurich  
Switzerland  
Tel. +41 44 709 01 01  
Fax +41 44 710 21 36  
[marflow@marflow.ch](mailto:marflow@marflow.ch)  
[www.marflow.ch](http://www.marflow.ch)

**MarFlow**<sup>®</sup>  
SWITZERLAND

# Nephrostomy Drainage Set with Pigtail Catheter

Made of Polyurethane  
Tecoflex® USA



# Nephrostomy Drainage Set with Pigtail Catheter

A convenient complete set for percutaneous placement of drainage catheter. Soft polyurethane material provides superior biocompatibility and drainage. Highly radiopaque for precise placement. Laying time up to 30 days

1 Set consists of:

- PIGTAIL – Drainage Catheter 35 cm, usable length 30 cm, plus straightener
- Puncture needle 1,3 mm (18 GA)
- Lunderquist Guidewire 80 cm J-Tip, fixed core, stainless steel, 0.038"
- Fascial Dilators 20 cm, ABS
- Urine Bag connector

## Set with PIGTAIL Catheter

Catheter hydrophilic coated

Art. No.	Catheter Ø (Fr/CH)	Fascial Dilators Ø (Fr/CH)
<b>PDS-P6</b>	6	5 – 6 – 7
<b>PDS-P7</b>	7	6 – 7 – 8
<b>PDS-P8</b>	8	7 – 8 – 9
<b>PDS-P9</b>	9	7 – 8 – 9 – 10
<b>PDS-P10</b>	10	6 – 8 – 10 – 12
<b>PDS-P12</b>	12	6 – 8 – 10 – 12
<b>PDS-P14</b>	14	6 – 8 – 12 – 14

## PIGTAIL Catheter only

<b>PCN-P6</b>	6
<b>PCN-P7</b>	7
<b>PCN-P8</b>	8
<b>PCN-P9</b>	9
<b>PCN-P10</b>	10
<b>PCN-P12</b>	12
<b>PCN-P14</b>	14



Marflow AG  
Soodstrasse 57  
8134 Adliswil/Zurich  
Switzerland  
Tel. +41 44 709 01 01  
Fax +41 44 710 21 36  
marflow@marflow.ch  
www.marflow.ch

May 2021

**MarFlow**<sup>®</sup>  
SWITZERLAND